


Drawings

US300PM Ultrasonic Flowmeter, Portable type US300PT Transducers for Portable type

Ultrasonic Flowmeter, Portable type
US300PM


Unit : mm
(approx. inch)



Weight: 3.9 kg (8.6 lb)


F01.EPS

Transducers for Portable type
US300PT-□B□



F02.EPS

Transducers for Portable type
US300PT-□C□
US300PT-□D□



F03.EPS

Transducers for portable type	X m (inch)	Y m (inch)	X+Y m (inch)	Weight kg (lb)
US300PT- □B□ (Medium type)	2.0 (78.7)	1.0 (39.4)	3.0 (118.1)	0.6 (1.32)
US300PT- □C□ (Large type)	2.0 (78.7)	2.4 (94.5)	4.4 (173.2)	1.2 (2.65)
US300PT- □D□ (Very large type)	5.0 (196.9)	7.0 (275.6)	12.0 (472.4)	2.2 (4.85)

T08-E.EPS