


Drawings


USPA301 or /WTG (for general temp.) USPA302 or /WTH (for high temp.) Wall Thickness probe

Unit : mm
(approx. inch)

Wall Thickness Probe for general temperature (−20° to +60°C / −4° to +140°F)
(model USPA301 or option /WTG)



Wall Thickness Probe for high temperature (0° to +200°C / 32° to +392°F)
(model USPA302 or option /WTH)



Length / Weight

Wall thickness probe	L m (inch)	Weight kg (lb)
USPA301 or /WTG	1.5 (59.1)	0.17 (0.37)
USPA302 or /WTH	1.2 (47.2)	0.19 (0.42)

T10.EPS