


CABLE LENGTH VARIES
 * PH97-07, 24 INCH CABLE LENGTH
 * PH97-21, 36 INCH CABLE LENGTH



PH87
 JUNCTION BOX
 TERMINAL STRIP

pH	CLEAR	# 5
REF	BLACK	# 3
TC	RED	# 1
TC	BLACK	# 2
SOLN GRND	GREEN	# 4

PH97 ELECTRODE
 SOLUTIONS GROUND, 0-14 pH GLASS
 PTFE REFERENCE JUNCTION
 Pt1000 TEMPERATURE SENSOR

UNITS:
 IN [MM]

YOKOGAWA

2 DART ROAD

NEWNAN, GEORGIA 30265