

## *FieldGuide INDEX* (Pressure)

<i>Topic</i>	<i>Transmitter</i>	<i>Protocol</i>	<i>Comm</i>	<i>Document #</i>	<i>Date</i>
Index				FGP-001A	January 2007
AL Option	<i>EJX</i>	Brain	BT200	FGP-101	February 2006
		Hart	275, 375	FGP-102	February 2006
Write Protect	<i>EJX</i>	Brain	BT200	FGP-111	February 2006
		Hart	275, 375	FGP-112	February 2006
Display	<i>EJX</i>	Brain	BT200	FGP-121	February 2006
		Hart	275, 375	FGP-122	February 2006
		Fieldbus	375-F	FGP-123	February 2006
High-Low Swap	<i>EJX</i>			FGP-130	February 2006
Response Time	<i>EJX</i>			FGP-140A	January 2007
Signal Characterization	<i>EJX</i>	Brain	BT200	FGP-151	March 2006
		Hart	375	FGP-152	March 2006
Calillary Fill Fluid	<i>EJX</i>	Brain	BT200	FGP-161	March 2006
Density Compensation		Hart	375	FGP-162	March 2006
Static Pressure	<i>EJX</i>	Hart	275, 375	FGP-182	January 2007
/EE Option	<i>EJA / EJX</i>	Fieldbus		FGP-501	January 2007